

## MAPKKide® Plus\*

### Specific Substrate for Anthrax Lethal Factor

*Bacillus anthracis* is regarded as a major biological warfare threat. The principal virulence factors are a γ-linked poly-D-glutamic acid (PGA) capsule and a three component exotoxin composed of protective antigen (PA), lethal factor (LF) and edema factor (EF). These proteins act in binary combinations. The complex of PA, the cell binding component, with the LF enzyme, is termed lethal toxin and can cause death.<sup>1,2</sup> PA and the enzymatic EF together cause skin edema.<sup>3</sup> Secreted PA is cleaved by membrane peptidases.<sup>4</sup> This allows the 63 kDa carboxy terminal fragment to oligomerize to a heptamer or higher.<sup>5</sup> Cleavage of PA is an essential step in exposing the binding sites for EF and LF. The complex enters the cell through endocytosis.<sup>6</sup> PA mediates the transfer of LF and EF to the cytoplasm where these enzymes recognize and alter their targets.

There are three forms of human anthrax: cutaneous, inhalation, and gastrointestinal. With the inhalational anthrax incident which resulted from the intentional release of anthrax spores in the bioterrorism attacks of 2001, early symptoms were similar to those for common illnesses.<sup>7</sup> Once symptoms are severe and diagnosis is possible, the levels of toxins are dangerously high. The incubation period can vary from hours to weeks depending on the dose received. Once intoxication occurs, anthrax bacteria can multiply rapidly in the blood and begin to secrete significant quantities of the toxin.<sup>8</sup> The inhalation form of *Bacillus anthracis* infection can kill quickly. While antibiotic treatment can clear the bacterium from the host, if diagnosis is delayed, the toxin, which is rapidly produced, may already be present in lethal amounts.

There is a critical need for a rapid, accurate, sensitive and simple assay to determine whether infection has occurred thereby allowing immediate treatment. Anthrax lethal factor (LF), an endopeptidase, is present in blood early in the infection.<sup>9</sup> Lethal factor is a zinc dependent metalloprotease which cleaves a specific bond in signaling proteins of the mitogen-activated protein kinase kinase family (MAPKK), destroying their ability to signal.<sup>10,11</sup> Of the seven different MAPKKs the amino terminus of six are cleaved by LF.<sup>12,13</sup> The crystal structure of LF complexed with the N-terminal portion of MAPKK-2 has been described.<sup>14</sup>

The use of peptidic substrates in plasma is problematic due to the presence of other proteases and the likelihood of nonspecific cleavage of the substrate. MAPKKide® Plus\*, Product #532 is a fluorogenic peptide substrate containing 7-amido-4-methylcoumarin (AMC) which is specific for LF, while resistant to nonspecific proteases. Significantly, this newly designed substrate is highly sensitive to LF and may be used to detect early infections. Two methods for the detection of LF in plasma have been developed using this fluorogenic substrate: an HPLC-based and a microplate assay.

Briefly, the LF is enriched by capture from plasma using an LF antibody-coated microtiter plate, and the captured LF is then exposed to the MAPKKide® Plus\* substrate. The amount of cleaved peptide substrate is determined by HPLC with fluorescence detection. Concentration of the LF using the antibody-coated plates allows for the detection of 5 pg LF/ml of neat plasma after 2 hours of incubation. Alternately the substrate may be added directly to diluted plasma and cleavage monitored by an increase in fluorescence as a function of time using a fluorescent microplate reader. The limit of detection by this simpler method is 1 ng LF/ml of plasma after 5 hours of digestion. Both methods can be confirmed by analysis of the reaction as a function of time.

MAPKKide® Plus\* (AMC) is supplied as a lyophilized powder, and a lot analysis detailing purity and a protocol describing the procedure for the specific detection of active infection of *Bacillus anthracis* in plasma accompanies each shipment.

\* Patent Pending

**These products are intended for research purposes only and are not for use in humans or as diagnostic agents. For further information, please contact List Biological Laboratories.**

©2016 LBL, Inc., Rev. 09/2017

Office: (408) 866-6363

Fax: (408) 866-6364

[www.listlabs.com](http://www.listlabs.com)

[info@listlabs.com](mailto:info@listlabs.com)



Related Products		
Product No.	Description	Size
530	MAPKKide® Peptide Substrate (o-Abz/Dnp) for <i>Bacillus anthracis</i> LF	200 nmoles
531	MAPKKide® Peptide Substrate (DABCYL/FITC) for <i>Bacillus anthracis</i> LF	200 nmoles
532	MAPKKide® Plus* (AMC) Specific Substrate for Anthrax Lethal Factor	100 nmoles
539	MAPKKide® Unquenched Calibration Peptide for Product #530	50 nmoles
169 A,B	Anthrax Lethal Factor (LF-A), Recombinant from <i>Bacillus anthracis</i> Native Sequence	0.1 mg, 1 mg
171 D,E	Anthrax Protective Antigen (PA), Recombinant from <i>Bacillus anthracis</i>	50 µg, 0.5 mg
172 B,C,D	Anthrax Lethal Factor (LF), Recombinant from <i>Bacillus anthracis</i>	70 µg, 1 mg, 10mg

\* Patent Pending

See how others have used List Labs' products on our citations page: <https://www.listlabs.com/citations>

## References

1. Leppla SH. **The bifactorial *Bacillus anthracis* lethal and oedema toxins.** In JE Alouf and JH Freer (ed.), The Comprehensive Sourcebook of Bacterial Protein Toxins. Second Edition, Academic Press, London. 1999; p. 243-263.
2. Mock M, Mignot T. Anthrax toxins and the host: a story of intimacy. *Cell Microbiol.* 2003; 5(1):15-23. [PMID:12542467](#)
3. Leppla SH. Anthrax toxin edema factor: a bacterial adenylate cyclase that increases cyclic AMP concentrations in eukaryotic cell. *Proc. Natl. Acad. Sci. U.S.A.* 1982; 79(10):3162-3166. [PMID:6285339](#)
4. Gordon VM, Klimpel KR, Arora N, Henderson MA, Leppla SH. Proteolytic activation of bacterial toxins by eukaryotic cells is performed by furin and by additional cellular proteases. *Infect Immun.* 1995; 63(1):82-87. [PMID:7806387](#)
5. Kintzer AF, Sterling HJ, Tang II, Williams ER, Krantz BA. Anthrax toxin receptor drives protective antigen oligomerization and stabilizes the heptameric and octameric oligomer by similar mechanism. *PLoS One.* 2010; 5(11):e13888. [PMID:21079738](#)
6. Singh Y, Klimpel KR, Goel S, Swain PK, Leppla SH. Oligomerization of anthrax toxin protective antigen and binding of lethal factor during endocytic uptake into mammalian cells. *Infect Immun.* 1999; 67(4):1853-1859. [PMID:10085027](#)
7. Jernigan JA, Stephens DS, Ashford DA, Omenaca C, Topiel MS, Galbraith M, Tapper M, Fisk TL, Zaki S, Popovic T, Meyer RF, Quinn CP, Harper SA, Fridkin SK, Sejvar JJ, Shepard CW, McDonnell M, Guarner J, Sheih WJ, Malecki JM, Gerberding JL, Hughes JM, Perkins BA. Bioterrorism-related inhalational anthrax: The first 10 cases reported in the United States. *Emerging Infect. Dis.* 2001; 7(6):933-943. [PMID:11747719](#)
8. Langer M, Duggan ES, Booth JL, Patel VI, Zander RA, Sllasi-Mansat R, Ramani V, Veres TZ, Prenzler F, Sewald K, Williams DM, Coggeshall KM, Awasthi S, Lupu F, Burian D, Ballard JD, Braun A, Metcalf JP. *Bacillus anthracis* lethal toxin reduces human alveolar epithelial barrier function. *Infect Immun.* 2012; 80(12):4374-4387. [PMID:23027535](#)

9. Boyer AE, Quinn CP, Hoffmaster AR, Kozel TR, Saile E, Marston CK, Percival A, Plikaytis BD, Woolfitt AR, Gallegos M, Sabourin P, McWilliams LG, Pirkle JL, Barr JR. Kinetics of lethal factor and poly-D-glutamic acid antigenemia during inhalation anthrax in rhesus macaques. *Infect Immun.* 2009; 77(8):3432-3441. [PMID:19506008](#)
10. Vitale G, Pellizzari R, Recchi C, Napolitani G, Mock M, Montecucco C. Anthrax lethal factor cleaves the N-terminus of MAPKKs and induces tyrosine/threonine phosphorylation of MAPKs in cultured macrophages. *Biochem. Biophys. Res. Commun.* 1998; 248(3):706-711. [PMID:9703991](#)
11. Pellizzari R, Guidi-Rontani C, Vitale G, Mock M, Montecucco C. Anthrax lethal factor cleaves MKK3 in macrophages and inhibits the LPS/IFNgamma-induced release of NO and TNF alpha. *FEBS Lett.* 1999; 462(1-2):199-204. [PMID:10580119](#)
12. Vitale G, Bernardi L, Napolitani G, Mock M, Montecucco C. Susceptibility of mitogen-activated protein kinase kinase family members to proteolysis by anthrax lethal factor. *Biochem. J.* 2000; 352:739-745. [PMID:11104681](#)
13. Pannifer AD, Wong TY, Schwarzenbacher R, Renatus M, Petosa C, Bienkowska J, Lacy DB, Collier RJ, Park S, Leppla SH, Hanna P, Liddington RC. Crystal structure of the anthrax lethal factor. *Nature*, 2001; 414(6860):229-233. [PMID:11700563](#)
14. Forino M, Johnson S, Wong TY, Rozanov DV, Savinov AY, Li W, Fattorusso R, Becattini B, Orry AJ, Jung D, Abagyan RA, Smith JW, Alibek K, Liddington RC, Pellecchia M. 2005, Efficient synthetic inhibitors of anthrax lethal factor. *Proc. Natl. Acad. Sci. U.S.A.* 2005; 102(27):499-9504. [PMID:15983377](#)

### Posters

1. Shine N, Eaton L, Crawford K. [Internally Quenched Fluorogenic Substrates for Anthrax Lethal Factor.](#) 2003.
2. Shine N, Suryadi K. [Fluorescence Substrates for Specific and Quantitative Detection of Anthrax Lethal Factor in Plasma.](#) 2016.