



LIST BIOLOGICAL LABORATORIES, INC.

540 DIVISION STREET ▪ CAMPBELL ▪ CALIFORNIA 95008-6906 ▪ USA
408-866-6363 ▪ 800-726-3213 ▪ FAX 408-866-6364 ▪ EMAIL info@listlabs.com
WEBSITE www.listlabs.com

Product #434
Date of Manufacture: April 2013
Retest Date: April 2018

CERTIFICATE OF ANALYSIS LIPOPOLYSACCHARIDE, ULTRA PURE from *Salmonella minnesota* R595 (Re) Lot #4345A1

Contents

Each vial contains 1 mg of lipopolysaccharide (LPS) lyophilized in water. This LPS is isolated from *Salmonella minnesota* R595 (Re) by a modification of the methods of Galanos *et al.*¹, and Manthey and Vogel.²

Reconstitution/Storage

LPS is dispersible in aqueous solvents at concentrations of 1 mg/ml. To achieve suspension in water, heating to about 50°C with intermittent vortexing or sonication is recommended.³ Allow ample time for dispersion to occur. The use of 0.5% triethylamine with brief vortexing aids in solubilization and dispersion. Triethylamine is very basic and may be neutralized with Tris HCl to avoid hydrolysis of the LPS fatty acid chains. It is recommended that this material be stored at 2-8°C prior to and following reconstitution.

Analysis

2-Keto-3-deoxyoctonate (KDO) ⁴	9.3%
Protein by Colloidal Gold Transfer blot ^{5,6}	≤0.01%
Nucleic acid ⁷	0.4%
Endotoxin activity by a kinetic chromogenic LAL assay.....	1.1 EU/ng

Handling

Good laboratory technique should be employed in the safe handling of any lipopolysaccharide or lipid A product. Wear appropriate laboratory attire including lab coat, gloves and safety glasses. Nitrile gloves are recommended when handling lyophilized material.

This product is pyrogenic. Avoid accidental autoinoculation by exercising extreme care when handling in conjunction with any injection device.

The product is intended for research purposes by qualified personnel. It is not intended for use in humans or as a diagnostic agent. List Biological Laboratories, Inc. is not liable for any damages resulting from the misuse or handling of this product.

FOR RESEARCH PURPOSES ONLY. NOT FOR HUMAN USE.

(continued)

Product #434, Lot #4345A1

References

1. Galanos, C., Lüderitz, O. and Westphal, O. (1969) *Eur. J. Biochem.* **9**, 245-249.
2. Manthey, C.L. and Vogel, S.N. (1994) *J. Endotoxin Research* **1**, 84-91.
3. Mukerjee, P., Kastowsky, M., Obst, S., Takayama, K. (1999) Lipopolysaccharide Preparations in Aqueous Media in *Endotoxin in Health and Disease*, Brade, H., Opal, S.M., Vogel, S.N., Morrison, D.C. eds., Marcel Dekker, Inc., New York, p. 223-224.
4. Cynkin, M.A. and Ashwell, G. (1960) *Nature* **186**, 155-156.
5. Rohringer, R., and Holden, D.W. (1985) *Anal. Biochem.* **144**, 118-127.
6. Danscher, G. (1981) *Histochemistry* **71**, 81-88.
7. Determined by absorption at 260 nm.

Production: VFK Date: 5/5/16

Management: NS Date: 5/5/16

QA/QC: TPP Date: 5/5/16

Made in U.S.A. 